

LiquiPure® Replacement Cartridges for Millipore® Laboratory Water Systems

LiquiPure® 1 Replacement Cartridges for Milli-Q® and Alpha-Q™ Laboratory Water Systems

The LiquiPure1 replacement cartridges for Milli-Q and Alpha-Q laboratory water systems provide ultra-high-purity water with very low TOC and extractable levels. The cartridge contains a variety of purification media to remove specific types of feed water contaminants. All wetted parts are constructed of all-natural polypropylene to minimize extractables, while a unique downflow design ensures optimum removal of organic and inorganic impurities. A pigmented polypropylene exterior prevents light from penetrating to minimize microbial growth. Cartridge contains macroreticular activated carbon, ultra-high-purity mixed bed ion exchange resin, and an organic scavenger.

ORDERING INFORMATION	Siemens Catalog Number	Millipore Equiv. Part Number
LiquiPure 1 Cartridge	LP1000	CPMQ004R1 CPMQ004D2
LiquiPure 1 Cartridge w/0.2 µm Final Filter	LP2000	CPMQK05R1 CPMQK05D2 CPMQK05RE
Final Filter 0.2 µm (1/pk)	LC134	MPGL04SK2

LiquiPure® CLC Replacement Cartridges for Millipore® Analyzer Feed Systems

The LiquiPure CLC replacement cartridges for Millipore AFS Analyzer Feed systems provide high-purity water to feed laboratory clinical chemistry analyzers. The cartridge contains ultra-high-purity mixed bed ion exchange resin to produce NCCLS Type 1 quality water. All wetted parts are constructed of all natural polypropylene to minimize extractables, while the unique design ensures optimum removal of inorganic impurities. A pigmented polypropylene exterior prevents light from penetrating to minimize microbial growth. Cartridge contains ultra-high-purity mixed bed ion exchange resin.

ORDERING INFORMATION	Siemens Catalog Number	Millipore Equiv. Part Number
LiquiPure CLC Cartridge	LP3000	CP4ALLRES
LiquiPure CLC Cartridge	LP3000R	CP4ALLRES
LiquiPure CLC Cartridge	LP3000R & LP4000 FCCFH11S2	CPKITAFS2
LiquiPure CLC Cartridge	LP3000R FCCFH11S2	CPAFS04R1

LiquiPure® RO Replacement Cartridges for Millipore Milli-RO® and AFS® Systems

These duplex RO pretreatment cartridges for Millipore Milli-RO and AFS Analyzer Feed systems includes granular activated carbon for chlorine removal, depth filtration for particulate removal and a food grade sequestering agent to protect the RO membrane from fouling and scale formation. Designed in a duplex configuration, the cartridge can be rotated and reconnected when one-half has been exhausted. All wetted parts are constructed of all natural polypropylene to minimize extractables, while the pigmented polypropylene exterior prevents light from penetrating to minimize microbial growth. Cartridge contains granular activated carbon, an anti-scaling agent and a prefilter.

ORDERING INFORMATION	Siemens Catalog Number	Millipore Equiv. Part Number
LiquiPure RO Cartridge	LP4000	CPR0P0402
LiquiPure RO Cartridge	LP4000R	CPR0P0402



Replacement Cartridges for Millipore® Laboratory Water Systems

ORDERING INFORMATION

Replacement Cartridges for Milli-Q® and Other Millipore Systems	Where Used	Millipore Equiv. Part Number	Siemens Catalog Number
High Capacity Carbon (4/pk)	3, 4, 5 Bowl Milli-Q or Milli-Q Plus Systems	CDFC01204	DICM1000-4
Mixed Bed (2/pk), Exchange Capacity, 540 Grains as CaCO ₃	3, 4, 5 Bowl Milli-Q or Milli-Q Plus Systems	CPMB01202	DIMM1200-2
Organic Scavenger (1/pk)	4 & 5 Bowl Milli-Q or Milli-Q Plus Systems	CDEX01201	DIOM1200-1
UF Membranes, 10k MW (1/pk)	4 & 5 Bowl Milli-Q or Milli-Q Plus Systems	CDUF01201	CDUF01201
Prefilter 12" 5 µm	Milli-RO® Systems	CDPRM1206	FCR001205
1 ea. DICM 1000-1 Carbon, 2 ea. DIMM 1200-2 Ion Exchange CWCF 1492 Capsule Final Filter (1ea.)	4-Bowl Life Science or 3-Bowl Standard	CFUF01204	LBCCM1202
1 ea. DICM 1000-1 Carbon, 2 ea. DIMM 1200-2 Ion Exchange, 1 ea. DIOM 1200-1 Organic Scavenger, 1 ea. CWCF 1492	4-Bowl Analytical System	CFOF01205	LBRCM1002
1 ea. DICM 1000-1 Carbon, 2 ea. DIMM 1200-2, 1 ea. CWCF 1492, 1 ea. 10" Polysulfone 0.2 µm Cartridge	5-Bowl Milli-Q	CFIF01205	LBCCM1202 FCWN010S2
1 ea. FCWN12152, DICM 1000-1, DIMM 1200-1, DIOM 1200-1, CWCF 1492	5-Bowl Milli-Q	CFOF012P5	FCWN010S2 DICM1000-1 DIMM1200-1 DIOM1200-1 CWCF 1492

Replacement Cartridges for Milli-Q® Ultrapure Water Systems	Where Used	Millipore Equiv. Part Number	Siemens Catalog Number
Organex Cartridge	Milli-Q Academic, Biocel, Synthesis, Element, Gradient	QTUM00EX	DIMOS1201
Ionic Cartridge	Milli-Q Academic, Biocel, Synthesis, Element, Gradient	QTUM00IX	DIMEX1201
Purification Pack (Elix, RO & Distillation Feed)	Milli-Q Academic, Biocel, Synthesis, Element, Gradient	QGARD00R1	PPMR01201
Purification Pack (DI Feed)	Milli-Q Academic, Biocel, Synthesis, Element, Gradient	QGARD00D2	PPMDI1201
Pyroguard 5000 MW Hollow Fiber UF Cartridge	Milli-Q Biocl, Synthesis	CDUFHF05K	CQUF500
VOC Cartridge	Milli-Q Academic, Biocel, Synthesis Element Gradient	QTUM000VX	DICMVOL01
Pretreatment Filter	Milli-Q Elix 3,5,10	PROG00002	PRCM1201
Pretreatment Filter	Milli-Q Elix 3,5,10	PROG00001	PRCMS001

Replacement Cartridge Kits for Millipore Super-Q® Systems	Where Used	Millipore Equiv. Part Number	Siemens Catalog Number
Activated Carbon (1/pk)	4-Bowl Super-Q (Housing #1)	CDFC02203	DICM2200-1
Ion Exchange (>2000 Grain Capacity) (1/pk)	4-Bowl Super-Q (Housing #2 & #3)	CPMB02202	DIMM2200-1
Organic Scavenger (1/pk)	4-Bowl Super-Q (Housing #4)	CDEX02201	DIOM2200-1
UF Cartridges	4-Bowl Super-Q (Housing #4 optional)	CDUFM2201	CDUF02201
0.22 µm Filter 20" (1/pk)	4-Bowl Super-Q (Housing #4 optional)	CVDI02TPE	FCDIO20S2

Replacement RO Membranes for Millipore Milli-RO® Plus and AFS® Systems	Where Used	Millipore Equiv. Part Number	Siemens Catalog Number
TFC RO Membrane	Milli-RO 4 & 15	CDRC012S0	ROGTM3512
TFC RO Membrane (w/out housing)	Milli-RO 6-10 and AFS Systems	CDRC012HG	LABROKIT
TFC RO Membrane (w/out housing)	Milli-RO 30/60/90	CDRCH12HG	LABROKIT
CA Membrane	Milli-RO 4 & 15	CDRO012S0 SL2M129V8	ROETM3512
TFC RO Membrane	Milli-RO 60/120	60 - CDRC025S1 120 [CDRC025S1 CDRC025SH	60 - CDRC025S1 120 [CDRC025S1 CDRC025SH

